

FITC Mouse IgG2a Isotype Antibody

Catalog Number	Vial Size
M20021-02B	50 µg
M20021-02D	200 µg
M20021-02G	25 tests
M20021-02H	100 tests



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
标准 高效 稳定 Precision Efficient Stable

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Clone	Isotype
C1.18.4	Mouse IgG2a

Description

The C1.18.4 immunoglobulin has unknown specificity. The isotype of this antibody is Mouse IgG2a,k. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

Product Information

Conjugation: FITC

Formulation: PBS pH 7.2, 0.09% NaN₃, 0.2% BSA

Storage: Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis.

For immunofluorescent staining, **For test sizes**, the suggested use of this reagent for immunofluorescent staining is 5 µl per million cells or 5 µl per 100 µl of whole blood.

For special application, the reagent has been pre-diluted for use at recommended Volume per test.

For µg sizes, use at concentrations comparable to those of the specific antibody of interest.

Application: Recommended Application: FC

References

- [1] Barclay, A.N., et al. 1997. The Leukocyte Antigen FactsBook, Academic Press.
- [2] Kishimoto, T.K., et al. 1990. P. Natl. Acad. Sci. USA 87:2244.
- [3] Tedder, T.F., et al. 1995. J. Exp. Med. 181:2259.

For Research Use Only.