

**IF555 Goat Anti-Rat IgG(H+L)**

Catalog Number	Vial Size
GR201G-37C	100 ug
GR201G-37E	500 ug



天津三箭生物技术股份有限公司  
Tianjin Sungene Biotech Co., Ltd.  
标准 高效 稳定 Precision Efficient Stable

**Market** | 400-621-0003  
marketing@sungenebiotech.com

**Support** | 022-66211636-8024  
techsupport@sungenebiotech.com

**Web** | www.sungenebiotech.com

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.  
This product is guaranteed up to one year from purchase.

**Purified Antibody Characterization**

Isotype	Reactivity
Polyclonal Goat IgG	Rabbit

**Description**

**Preparation:** Purified from goat antiserum by rat IgG affinity chromatography.

**Purity:** >90% by SDS-PAGE

**Specificity:** The antibody reacts with whole molecule Rat IgG. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.

**Product Information**

**Conjugation:** IF555

**Formulation:** PBS pH 7.2, 0.09% NaN<sub>3</sub>, 0.2% BSA

**Concentration:** 0.5 mg/ml

**Storage:** Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

**Application:** IF(Immunofluorescent)

**Usage:** IF:1/10-1/1000

**The appropriate dilutions must be determined for individual use.**

**For Research Use Only.**